





Product

phenofAST Product



Disposable cassette for testing multiple antibiotics for each patient



Each device process 4 cassettes. Labs can stack multiple devices to increase throughput and save space



PhenofAST provides automatic antibiotic susceptibility report

Providers survey: Market values are
accuracy, time-to-result, cost-effectiveness
and broad performance.

Gatekeepers	Frontline Clinicians (EM, ICU, ID)	Microbiologist Laboratory Directors	Operational Laboratory Directors
Priorities for AST	Accuracy Time-to-results Clinically relevant antibiotics	Accuracy Time-to-results Wide range test menu	Operational efficiency Cost effectiveness



FOR FURTHER DETAILS

Eliane Markowitz PhD

+972-54-492-0616 | eliane.m@phenofast.com

Nadav Ben Dov PhD

+972-50-590-5060 | nadav.bd@phenofast.com